






# Heart Assist Lab (BME Cornell)

## CVS Simulator Manual (Version 3)

Hiral Kothari, James F. Antaki

This is the 3<sup>rd</sup> version of the cardiovascular (CVS) simulator, developed by the Heart Assist Lab at Cornell University. It is aimed to simulate the hemodynamics of human cardiovascular and provide a non-invasive way to investigate the intervention of heart disease, as well as the intervention of rotary ventricular assist device for the patients with end-stage heart failure. This model includes bi-ventricles, systemic circulation, pulmonary circulation, the regulation of normal baro-reflex<sup>1</sup>, biventricular interaction<sup>2-4</sup> and model of rotary left ventricular assist device (VAD)<sup>5</sup>. It enables the user to simulate different scenarios of heart diseases and VAD therapy by changing the set up in each tab (Please see instruction 1-4). We also provide five examples to facilitate your beginning trial of this CVS simulator (Please see the Instructions at the end).

When you finish installation, you'll find a short-cut of the software on your desktop. (Note: make sure you put the installation folder under the C:\ path, since the program links an excel file in the folder and set the path as "C:\CVS Installer\data.xls").

Double click the short-cut; the program begins to run, when you click  on the up-left corner of the interface. You could click the  symbol on the up-left  corner or on the interface to end the simulation. For all the control panels and waveform display, you could adjust the scale of the axis by clicking the maximum or minimum values to input other digits.

This has the following features:

- There are 3 tabs for easy access of the controls and the indicators: User Guide, Controls & Monitor, Dilated HF and Septum Monitor. All the controls are colored yellow, and all the indicators are in pink.
- A data file with for easier understanding.



# Heart Assist Lab (BME Cornell)

**User Guide (Figure 1):** Brief instructions to get you started right away. This tab will help you navigate through the software and provide important guidelines to avoid errors.





| Reference Set Up   | Before You Start  |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
|--|-------------------|-----------------|--------------------------------|-----|--------------------------------|-------|----------------------|-----|--|---|---|------|---|-------|------------------------------|------|----------------|--|---------------|-------|---------------|----|---------------|------|---------------|-----|-------------|-----|-------------|-----|-----------|-------|-----------|------|---|
| <p>For a normal cardiovascular system :</p> <p>Select the "All Baro" option at the Baroreflex Control and <b>turn off</b> the Pump.</p> <p>Set the parameters as :</p> <table border="1"><thead><tr><th>For Heart Disease</th><th>Normal Baseline</th></tr></thead><tbody><tr><td>LV Elastance (Elvmax: mmHg/ml)</td><td>2.4</td></tr><tr><td>RV Elastance (Ervmax: mmHg/ml)</td><td>1.412</td></tr><tr><td>Heart Rate (HR: bpm)</td><td>103</td></tr><tr><td>Stenosis (AV, MV, TV, PV: mmHg<sub>soml</sub>-1)</td><td>0</td></tr><tr><td>LV afterload (SVR: mmHg<sub>soml</sub>-1)</td><td>0.06</td></tr><tr><td>RV afterload (PVR: mmHg<sub>soml</sub>-1)</td><td>0.023</td></tr><tr><td>Total Blood Volume (TBV: ml)</td><td>5300</td></tr></tbody></table><br><table border="1"><thead><tr><th>For Dilated HF</th><th></th></tr></thead><tbody><tr><td>Csv (ml/mmHg)</td><td>61.11</td></tr><tr><td>Cev (ml/mmHg)</td><td>50</td></tr><tr><td>Cpa (ml/mmHg)</td><td>0.76</td></tr><tr><td>Cpp (ml/mmHg)</td><td>5.8</td></tr><tr><td>Polv (mmHg)</td><td>1.5</td></tr><tr><td>Porv (mmHg)</td><td>1.5</td></tr><tr><td>Vulv (ml)</td><td>16.77</td></tr><tr><td>Vurv (ml)</td><td>46.8</td></tr></tbody></table> | For Heart Disease | Normal Baseline | LV Elastance (Elvmax: mmHg/ml) | 2.4 | RV Elastance (Ervmax: mmHg/ml) | 1.412 | Heart Rate (HR: bpm) | 103 | Stenosis (AV, MV, TV, PV: mmHg <sub>soml</sub> -1) | 0 | LV afterload (SVR: mmHg <sub>soml</sub> -1) | 0.06 | RV afterload (PVR: mmHg <sub>soml</sub> -1) | 0.023 | Total Blood Volume (TBV: ml) | 5300 | For Dilated HF |  | Csv (ml/mmHg) | 61.11 | Cev (ml/mmHg) | 50 | Cpa (ml/mmHg) | 0.76 | Cpp (ml/mmHg) | 5.8 | Polv (mmHg) | 1.5 | Porv (mmHg) | 1.5 | Vulv (ml) | 16.77 | Vurv (ml) | 46.8 | <p>Make sure you put the installation folder under the C:\ path, since the program links a excel file in the folder and set the path as "C:\CVS_V2_Installer\data.xls.</p> <p>If an error pops up after this, then go to the <b>Data Record</b> section and click on the <b>destination path</b>. Manually insert the .xls file named "data".</p> <p>Click on the  symbol on the upper left corner to start the simulation and on symbols  at the up-left corner or  on the interface to stop the simulation.</p> <p>To adjust the scale of the axis on the controls and waveform charts, click the maximum or minimum values to input other valu .</p> <p>For accurate Ramp Speed visuals, select the Ramp speed first. Shifting from Manual/Constant Speed to Ramp creates undesired results.</p> <p style="text-align: center;"><b>Please refer to the User Manual pdf for detailed description and instructions</b></p> |
| For Heart Disease  | Normal Baseline   |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| LV Elastance (Elvmax: mmHg/ml)   | 2.4               |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| RV Elastance (Ervmax: mmHg/ml)   | 1.412             |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Heart Rate (HR: bpm)   | 103               |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Stenosis (AV, MV, TV, PV: mmHg <sub>soml</sub> -1)   | 0                 |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| LV afterload (SVR: mmHg <sub>soml</sub> -1)  | 0.06              |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| RV afterload (PVR: mmHg <sub>soml</sub> -1)  | 0.023             |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Total Blood Volume (TBV: ml)   | 5300              |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| For Dilated HF   |                   |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Csv (ml/mmHg)  | 61.11             |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Cev (ml/mmHg)  | 50                |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Cpa (ml/mmHg)  | 0.76              |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Cpp (ml/mmHg)  | 5.8               |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Polv (mmHg)  | 1.5               |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Porv (mmHg)  | 1.5               |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Vulv (ml)  | 16.77             |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| Vurv (ml)  | 46.8              |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |
| <p>With the baroreflex regulation, the changed values of Heart rate, LV Elastance and RV Elastance can be seen in the baroreflex monitor.</p> <p>Only the combination of "linear" LV/RV ESPVR hypothesis in Dilated HF works under "All Baro"; all the other combinations have to work under either "No Baro" or "Non-heart" setting.</p>  |                   |                 |                                |     |                                |       |                      |     |  |   |   |      |   |       |                              |      |                |  |               |       |               |    |               |      |               |     |             |     |             |     |           |       |           |      |   |

Figure 1: User Guide Tab

**Controls and Monitors 1 (Figure 2):** Settings for general Heart Disease, LVAD controls, Baroreflex Control and Cardiac Monitor.

1. **Heart Diseases:** This can be created by manipulating the following-

- Left Ventricle (LV)/ Right Ventricle (RV) contractility based on the hypothesis of the linear end-systolic pressure and volume relationship (ESPVR) (mmHg/mL)  
*Note: "Dilated HF" default settings for LV/RV ESPVR Hypothesis are linear.*
- LV/RV afterload (e.g. systemic/pulmonary hypertension) (mmHg\*s\*mL<sup>1</sup>)
- Stenosis of Aortic Valve (AV), Mitral Valve (MV), Tricuspid Valve (TV) & Pulmonary Valve (PV) (mmHg\*s\*mL<sup>-1</sup>)
- Regurgitation of Aortic Valve and Mitral Valve
- Heart rate (HR) (bpm)
- The total blood volume (hypervolemia/hypovolemia) (mL)




## Heart Assist Lab (BME Cornell)

2. **LVAD (Rotary VAD) Controls:** You could turn on the LVAD in heart failure condition to see its intervention on cardiovascular hemodynamics. There are three options for the speed setup:

- Ramp Speed: It begins at 2k and goes up to 16k. You could get the instantaneous Rotary VAD speed from the meter of “*Pump Speed*”.
- Manual Speed: you can set up the speed from dragging the slider from 2k to 12k.
- Pulsatile Speed: The rotating speed of the VAD uses the sine function to create pulses. **Switch on the Pulsatile Speed button after switching to constant/manual speed.** This feature does not work with Ramp Speed

All these controls are grouped under the name of Rotary VAD.

3. **Data Record:** This session provides you an option to record the data. To create an excel file on your computer click the symbol  and select the directory named “data” to save the recorded data. If you do not select the directory correctly, when you stop the program, there will be an error showing up.

4. **Baroreflex Control:** There are three control options to maintain 7 parameters of the model to maintain arterial pressure. The barometer monitor shows the instantaneous values of the 7 parameters running in the program.

- “No Baro”: no baro-reflex regulation – “baro-reflex” control is off.
- “All Baro”: regulates Heart Rate, Contractility of LV, RV free wall, Peripheral Systemic Resistance, including both splanchnic  $R_{sp}$  and extra-splanchnic vessels ( $R_{ep}$ ) and Systemic Venous Unstressed Volume, including both splanchnic  $V_{usv}$  and extra-splanchnic vessels  $V_{uev}$ . (Note: only the combination of “linear” LV/RV ESPVR hypothesis works under “All Baro”; all the other combinations must work under either “No Baro” or “Non-heart” setting)
- “Non-heart”: regulates  $R_{sp}$ ,  $R_{ep}$ ,  $V_{usv}$  and  $V_{uev}$  (default settings)



# Heart Assist Lab (BME Cornell)




Figure 2: Controls and Monitor Tab

**Dilated HF and Septum Monitor (Figure 3):** Dilated Heart Failure controls, Mean value for certain parameters and instantaneous monitor for septum <sup>4</sup>.

1. **Dilated Heart Failure (HF):** You could set the four heart failure conditions listed below or manipulate the parameters according to your needs.
  - **CVS Setting:** Here you can change the LV/RV ESPVR hypothesis - “linear” refers to traditional linear ESPVR settings (see the dotted line in Figure 4 labeled with normal); all the other “nonlinear” conveys the concept that when the cardiac myosin overstretches, the force produced by them will decrease (see the dotted line in Figure 3 labeled with depressed ESPVR).



# Heart Assist Lab (BME Cornell)

The following are the equations to model LV and RV; Table 1 lists the settings for each of the nonlinear ESPVR hypothesis.

$$P_{\max,lv}(t) = e(t) \cdot P_{lv,es} + [1 - e(t)] \cdot P_{ed,lv} \quad 0 \leq e(t) \leq 1 \quad 1$$

$$e(t) = \begin{cases} \sin^2(\pi \cdot T(t) \cdot u / T_{sys}) & 0 \leq u \leq T_{sys} / T \\ 0 & T_{sys} / T \leq u \leq 0 \end{cases} \quad 2$$

$$P_{es,lv} = \begin{cases} 0 & V_{lv} \leq V_{u,lv} \\ E_{+,lv}(V_{lv} - V_{u,lv}) & V_{u,lv} < V_{lv} \leq (E_{+,lv}V_{u,lv} + V_{b,lv}) / (E_{+,lv} - E_{-,lv}) \\ E_{-,lv}V_{lv} + V_{b,lv} & V_{lv} > (E_{+,lv} \cdot V_{u,lv} + V_{b,lv}) / (E_{+,lv} - E_{-,lv}) \end{cases} \quad 3$$

$$P_{ed,lv} = P_{0,lv} \cdot (e^{k_{E,lv}(V_{lv} - V_{u,lv})} - 1) \quad 4$$

$$P_{lv} = P_{\max,lv} - R_{lv} \cdot Q_{o,l} \quad 5$$

$$R_{lv} = k_{R,lv} \cdot P_{\max,lv} \quad 6$$

**Table 1: Settings for each of the nonlinear ESPVR hypothesis**

|                                  | $E_{+,lv}$ | $V_{u,lv}$ | $E_{-,lv}$ | $V_{b,lv}$ |
|----------------------------------|------------|------------|------------|------------|
| Heart failure (nonlinear ESPVR1) | 1.2        | 25         | -0.6       | 240        |
| Heart failure (nonlinear ESPVR2) | 1          | 25         | -0.6       | 215        |
| Heart failure (nonlinear ESPVR3) | 0.8        | 25         | -0.6       | 190        |
| RV failure (nonlinear ESPVR1)    | 9/8        | 40.8       | -0.2       | 117.36     |
| RV failure (nonlinear ESPVR2)    | 0.9        | 40.8       | -0.1       | 73.28      |
| RV failure (nonlinear ESPVR3)    | 13/8       | 40         | -0.4       | 97         |
| RV failure (nonlinear ESPVR4)    | 1.6        | 40.8       | -0.2       | 204.72     |



# Heart Assist Lab (BME Cornell)

LV/RV input provides the options to adjust the LV/RV diastolic functions.

You could also change the parameters in the systemic/pulmonary circulation to simulate different vascular diseases.

- **Baroreflex Setting:** This offers the opportunity to adjust the parameters in the baroreflex system, including “unstressed volume”, “Vusv Baro Parameter” and “Vuev Baro Parameter”. More detailed definitions for each parameter are in Ursino’s study <sup>1</sup>.






Figure 3: Dilated HF and Septum Monitor



# Heart Assist Lab (BME Cornell)



(a)



(b)

Figure 4: Pressure-volume relationships for normal and four heart failure conditions: (a) left ventricle; (b) right ventricle.





# Heart Assist Lab (BME Cornell)

**Reference Set-Up:** Here are the few suggested settings for your trial of simulation: -

- 1 For one normal cardiovascular system, we select “All Baro” baro-reflex control and turn off the Rotary VAD; set the parameters as the ones listed in the Table 2 (normal baseline condition). With the Baroreflex regulation, the real value of HR, LV free wall elastance and RV free wall elastance regulated in the program are displayed in the “Baroreflex monitoring” panel.
- 2 For another normal cardiovascular system, we select “Non-heart Baro” control and turn off the Rotary VAD; set the parameters as the one listed in the Table 4 (the default setting for this program is to simulate these conditions and the corresponding PV loop, the hemodynamic outputs are detailed in Table 7.)
- 3 For 4 heart failure conditions, we select “Non-heart Baro” control in the “Baroreflex” tab; set all the parameters listed in Table 6 to the indicated values in each condition.
- 4 Table 7 demonstrates the corresponding hemodynamic output under each of the four conditions.
  - LV failure: LV ESPVR hypothesis “heart failure (nonlinear ESPVR1)”; RV ESPVR hypothesis “linear” or “RV normal (nonlinear ESPVR)”.
  - Bi-ventricular failure: LV ESPVR hypothesis “heart failure (nonlinear ESPVR1)”; RV ESPVR hypothesis “RV failure (nonlinear ESPVR1)”.
  - LV failure with PVR+: LV ESPVR hypothesis “heart failure (nonlinear ESPVR1)”; RV ESPVR hypothesis “linear” or “RV normal (nonlinear ESPVR)”.
  - Bi-ventricular failure with PVR+: LV ESPVR hypothesis “heart failure (nonlinear ESPVR1)”; RV ESPVR hypothesis “RV failure (nonlinear ESPVR1)”.





# Heart Assist Lab (BME Cornell)

**Table 2: Reference set up for normal conditions with linear ESPVR hypothesis (default).**

| <b>“Heart Disease” Tab</b>                          | <b>Normal baseline</b> |
|---|------------------------|
| LV Elastance (Elvmax: mmHg/ml)                      | <b>2.4</b>             |
| RV Elastance (Ervmax: mmHg/ml)                      | 1.412                  |
| Heart Rate (HR: bpm)                                | <b>103</b>             |
| Stenosis (AV, MV, TV, PV: mmHg•s•ml <sup>-1</sup> ) | 0                      |
| LV afterload (SVR: mmHg•s•ml <sup>-1</sup> )        | 0.06                   |
| RV afterload (PVR: mmHg•s•ml <sup>-1</sup> )        | 0.023                  |
| Total Blood Volume (TBV: ml)                        | 5300                   |
| <b>“Dilated HF” Tab</b>                             | <b>Normal baseline</b> |
| Csv (ml/mmHg)                                       | 61.11                  |
| Cev (ml/mmHg)                                       | 50                     |
| Cpa (ml/mmHg)                                       | 0.76                   |
| Cpv (ml/mmHg)                                       | 5.8                    |
| Polv (mmHg)   | <b>1.5</b>             |
| Porv (mmHg)   | 1.5                    |
| Vulv (ml)   | <b>16.77</b>           |
| Vurv (ml)   | 46.8                   |



# Heart Assist Lab (BME Cornell)

**Table 3: The normal Parameters for the Vascular System**

| Parameter Group Name                    | Symbol     | Units                   | Values         |
|---|------------|-------------------------|----------------|
| <b><i>Compliance</i></b>                |            |                         |                |
| Systemic arterial circulation           | $C_{sa}$   | ml/mmHg                 | 0.28           |
| Splanchnic peripheral circulation       | $C_{sp}$   | ml/mmHg                 | 2.05           |
| Extra splanchnic peripheral circulation | $C_{ep}$   | ml/mmHg                 | 1.67           |
| Splanchnic venous circulation           | $C_{sv}$   | ml/mmHg                 | 61.11          |
| Extra splanchnic venous circulation     | $C_{ev}$   | ml/mmHg                 | 50.0           |
| Pulmonary arterial circulation          | $C_{pa}$   | ml/mmHg                 | 0.76           |
| Pulmonary peripheral circulation        | $C_{pp}$   | ml/mmHg                 | 5.80           |
| Pulmonary venous circulation            | $C_{pv}$   | ml/mmHg                 | 25.37          |
| <b><i>Unstressed Volume</i></b>         |            |                         |                |
| Systemic arterial circulation           | $V_{u,sa}$ | ml                      | 0              |
| Splanchnic peripheral circulation       | $V_{u,sp}$ | ml                      | 274.4          |
| Extra splanchnic peripheral circulation | $V_{u,ep}$ | ml                      | 336.6          |
| Splanchnic venous circulation           | $V_{u,sv}$ | ml                      | 1121           |
| Extra splanchnic venous circulation     | $V_{u,ev}$ | ml                      | 1375           |
| Pulmonary arterial circulation          | $V_{u,pa}$ | ml                      | 0              |
| Pulmonary peripheral circulation        | $V_{u,pp}$ | ml                      | 123            |
| Pulmonary venous circulation            | $V_{u,pv}$ | ml                      | 120            |
| <b><i>Hydraulic Resistance</i></b>      |            |                         |                |
| Systemic arterial circulation           | $R_{sa}$   | mmHg*s*ml <sup>-1</sup> | 0.06           |
| Splanchnic peripheral circulation       | $R_{sp}$   | mmHg*s*ml <sup>-1</sup> | 3.307          |
| Extra splanchnic peripheral circulation | $R_{ep}$   | mmHg*s*ml <sup>-1</sup> | 1.407          |
| Splanchnic venous circulation           | $R_{sv}$   | mmHg*s*ml <sup>-1</sup> | 0.038          |
| Extra splanchnic venous circulation     | $R_{ev}$   | mmHg*s*ml <sup>-1</sup> | 0.016          |
| Pulmonary arterial circulation          | $R_{pa}$   | mmHg*s*ml <sup>-1</sup> | 0.023          |
| Pulmonary peripheral circulation        | $R_{pp}$   | mmHg*s*ml <sup>-1</sup> | 0.0894         |
| Pulmonary venous circulation            | $R_{pv}$   | mmHg*s*ml <sup>-1</sup> | 0.0056         |
| <b><i>Inertance</i></b>                 |            |                         |                |
| Systemic arterial circulation           | $L_{sa}$   | mmHg*ml*s <sup>-2</sup> | $0.22*10^{-3}$ |
| Pulmonary arterial circulation          | $L_{pa}$   | mmHg*ml*s <sup>-2</sup> | $0.18*10^{-3}$ |
| <b><i>Total Blood Volume</i></b>        | <b>TBV</b> | <b>ml</b>               | <b>5300</b>    |



# Heart Assist Lab (BME Cornell)

**Table 4: The normal Parameters for left and right heart**

| Parameter Group Name                        | Symbol     | Units                   | Values                |
|---|------------|-------------------------|-----------------------|
| <b><i>Left Heart</i></b>                    |            |                         |                       |
| Left atrial compliance                      | $C_{la}$   | ml/mmHg                 | 19.23                 |
| Left atrial unstressed blood volume         | $V_{u,la}$ | ml                      | 16.7                  |
| Left atrial resistance                      | $R_{la}$   | mmHg*s*ml <sup>-1</sup> | 2.5*10 <sup>-3</sup>  |
| Scaling factor for LV EDPVR                 | $P_{0,lv}$ | mmHg                    | 1.5                   |
| Exponent for LV EDPVR                       | $K_{E,lv}$ | ml <sup>-1</sup>        | 0.017                 |
| LV unstressed blood volume                  | $V_{u,lv}$ | ml                      | 16.7                  |
| Ascending slope for LV ESPVR                | $E_{+,lv}$ | mmHg/ml                 | 2.9                   |
| Descending slope for LV ESPVR               | $E_{-,lv}$ | mmHg/ml                 | --                    |
| Descending LV ESPVR pressure axis intercept | $V_{b,lv}$ | ml                      | --                    |
| Scaling factor for LV viscous resistance    | $K_{R,lv}$ | s/ml                    | 3.75*10 <sup>-4</sup> |
| <b><i>Right Heart</i></b>                   |            |                         |                       |
| Right atrial compliance                     | $C_{ra}$   | ml/mmHg                 | 31.25                 |
| Right atrial unstressed blood volume        | $V_{u,ra}$ | ml                      | 25                    |
| Right atrial resistance                     | $R_{ra}$   | mmHg*s*ml <sup>-1</sup> | 2.5*10 <sup>-3</sup>  |
| Scaling factor for RV EDPVR                 | $P_{0,rv}$ | mmHg                    | 1.5                   |
| Exponent for RV EDPVR                       | $K_{E,rv}$ | ml <sup>-1</sup>        | 0.016                 |
| RV unstressed blood volume                  | $V_{u,rv}$ | ml                      | 40.8                  |
| Ascending slope for RV ESPVR                | $E_{+,rv}$ | mmHg/ml                 | 1.75                  |
| Descending slope for RV ESPVR               | $E_{-,rv}$ | mmHg/ml                 | --                    |
| Descending RV ESPVR pressure axis intercept | $V_{b,rv}$ | ml                      | --                    |
| Scaling factor for RV viscous resistance    | $R_{sa}$   | s/ml                    | 1.4*10 <sup>-3</sup>  |

\* LV: left ventricle; RV: right ventricle; ESPVR: End-systolic pressure and volume relationship; EDPVR: End-diastolic pressure and volume relationship



# Heart Assist Lab (BME Cornell)

**Table 5: The parameters for baroreflex control**

| Parameter name/description   | Symbol         | Units                                  | Values |
|--|----------------|--|--------|
| <b><i>Carotid sinus afferent pathway</i></b>                       |                |  |        |
| Intracranial pressure at the central point of sigmoidal function   | $P_n$          | mmHg                                   | 92     |
| Lower saturation of frequency discharge in afferent fibers         | $f_{\min}$     | spikes/s                               | 2.52   |
| Upper saturation of frequency discharge in afferent fibers         | $f_{\max}$     | spikes/s                               | 47.78  |
| Relate to the slope of the static function at the central point    | $k_a$          | mmHg                                   | 11.758 |
| Time constant for real zero in the linear dynamic block            | $\tau_z$       | s                                      | 6.37   |
| <b><i>Sympathetic efferent pathway</i></b>                         |                |  |        |
| Frequency constant of spikes in the efferent sympathetic nerves    | $f_{es,oo}$    | spikes/s                               | 2.10   |
| Frequency constant of spikes in the efferent sympathetic nerves    | $f_{es,0}$     | spikes/s                               | 16.11  |
| Minimum sympathetic stimulation                                    | $f_{es,\min}$  | spikes/s                               | 2.66   |
| Constant in the exponential frequency trend                        | $k_{es}$       | spikes/s                               | 0.0675 |
| <b><i>Vagal efferent pathway</i></b>                               |                |  |        |
| Frequency constant of spikes in the efferent vagal nerves          | $f_{ev,oo}$    | spikes/s                               | 6.3    |
| Frequency constant of spikes in the efferent vagal nerves          | $f_{ev,0}$     | spikes/s                               | 3.2    |
| Constant in the exponential frequency trend                        | $k_{ev}$       | spikes/s                               | 7.06   |
| Central value for the frequency of spikes in the afferent fibers   | $f_{cs,0}$     | spikes/s                               | 25     |
| <b><i>Effectors</i></b>  |                |  |        |
| Mechanism strength for splanchnic peripheral resistance            | $G_{R,sp}$     | mmHg•ml <sup>-1</sup> •v <sup>-1</sup> | 0.0695 |
| Mechanism strength for Extrasplanchnic peripheral resistance       | $G_{R,ep}$     | mmHg•ml <sup>-1</sup> •v <sup>-1</sup> | 0.53   |
| Mechanism strength for unstressed splanchnic venous volume         | $G_{Vu,sv}$    | ml•v <sup>-1</sup>                     | -265.4 |
| Mechanism strength for unstressed Extrasplanchnic venous volume    | $G_{Vu,ev}$    | ml•v <sup>-1</sup>                     | -132.5 |
| Time constant for splanchnic peripheral resistance                 | $\tau_{R,sp}$  | s                                      | 6      |
| Time constant for Extrasplanchnic peripheral resistance            | $\tau_{R,ep}$  | s                                      | 6      |
| Time constant for unstressed splanchnic venous volume              | $\tau_{Vu,sv}$ | s                                      | 20     |
| Time constant for unstressed Extrasplanchnic venous volume         | $\tau_{Vu,ev}$ | s                                      | 20     |
| Time delay for splanchnic peripheral resistance                    | $D_{R,sp}$     | s                                      | 6      |
| Time delay for Extrasplanchnic peripheral resistance               | $D_{R,ep}$     | s                                      | 2      |
| Time delay for unstressed splanchnic venous volume                 | $D_{Vu,sv}$    | s                                      | 5      |
| Time delay for unstressed Extrasplanchnic venous volume            | $D_{Vu,ev}$    | s                                      | 5      |
| Splanchnic peripheral resistance in absence of innervation         | $R_{sp,0}$     | mmHg•s•ml <sup>-1</sup>                | 2.49   |
| Extrasplanchnic peripheral resistance in absence of innervation    | $R_{ep,0}$     | mmHg•s•ml <sup>-1</sup>                | 0.78   |
| Unstressed splanchnic venous volume in absence of innervation      | $V_{usv,0}$    | ml                                     | 1435.4 |
| Unstressed Extrasplanchnic venous volume in absence of innervation | $V_{usev,0}$   | ml                                     | 1537   |



# Heart Assist Lab (BME Cornell)

**Table 6: Parameter values to simulate four pathological conditions**

|                               | <i>LVF</i> | <i>Bi-F</i> | <i>LVF, PVR+</i> | <i>Bi-F, PVR+</i> |
|-------------------------------|------------|-------------|------------------|-------------------|
| <b><i>Left Ventricle</i></b>  |            |             |                  |                   |
| $K_{E,lv}$                    | 0.0136     | 0.0144      | 0.0121           | 0.0109            |
| $V_{u,lv}$                    | 25         | 25          | 25               | 25                |
| $E_{+,lv}$                    | 1.2        | 1.2         | 1.2              | 1.2               |
| $E_{-,lv}$                    | -0.6       | -0.6        | -0.6             | -0.6              |
| $V_{b,lv}$                    | 240        | 240         | 240              | 240               |
| <b><i>Right Ventricle</i></b> |            |             |                  |                   |
| $K_{E,rv}$                    | 0.016      | 0.0115      | 0.016            | 0.0115            |
| $E_{+,rv}$                    | 1.6        | 0.8         | 1.6              | 0.8               |
| $E_{-,rv}$                    | -0.2       | -0.2        | -0.2             | -0.2              |
| $V_{b,rv}$                    | 204.72     | 117.36      | 204.72           | 117.36            |
| <b><i>Vascular System</i></b> |            |             |                  |                   |
| $R_{pa}$                      | 0.023      | 0.092       | 0.092            | 0.092             |
| $R_{pp}$                      | 0.0894     | 0.0894      | 0.44             | 0.44              |
| $V_{u,sp}$                    | 274.4      | 94.4        | 274.4            | 94.4              |
| $V_{u,ep}$                    | 336.6      | 66.6        | 336.6            | 66.6              |
| $V_{u,sv}$                    | 1435.4     | 735.4       | 1435.4           | 735.4             |
| $V_{u,ev}$                    | 1537       | 467         | 1537             | 467               |
| TBV                           | 5700       | 5700        | 5700             | 5700              |

Note that model parameter value not shown in this table retained normal values, as specified in Appendix 4. See Appendix 4 for definitions of abbreviations.

LVF: left ventricular failure; Bi-F: biventricular failure, including LVF and right ventricular failure (RVF); PVR+: high pulmonary vascular resistance.



## Heart Assist Lab (BME Cornell)

**Table 7: Simulation variables at normal and four pathological heart failure conditions.**

| <i>Parameters</i>             | <i>Units</i> | <i>Normal</i> | <i>LVF</i> | <i>Bi-F</i> | <i>LVF,PVR+</i> | <i>PVR+</i> |
|-------------------------------|--------------|---------------|------------|-------------|-----------------|-------------|
| Cardiac output                | L/min        | 5.8           | 4.5        | 3.2         | 3.3             | 2.4         |
| Stroke volume                 | mL           | 68.0          | 41.1       | 37.1        | 39.2            | 28.6        |
| Left atrial pressure (mean)   | mmHg         | 8.7           | 26.0       | 28.0        | 19.1            | 18.0        |
| Right atrial pressure (mean)  | mmHg         | 4.4           | 3.5        | 17.1        | 5.1             | 20.7        |
| Arterial pressure             |              |               |            |             |                 |             |
| Systolic                      | mmHg         | 133.0         | 95.2       | 98.4        | 93.0            | 87.6        |
| Diastolic                     | mmHg         | 93.0          | 69.4       | 74.8        | 68.3            | 68.0        |
| Mean                          | mmHg         | 106.0         | 77.3       | 82.1        | 75.9            | 73.8        |
| Pulmonary arterial pressure   |              |               |            |             |                 |             |
| Systolic                      | mmHg         | 28.0          | 37.9       | 49.8        | 61.9            | 49.4        |
| Diastolic                     | mmHg         | 15.0          | 29.3       | 31.7        | 42.3            | 35.0        |
| Mean                          | mmHg         | 20.1          | 32.9       | 37.8        | 48.9            | 39.8        |
| LV end-diastolic pressure     | mmHg         | 10.0          | 25.6       | 27.6        | 18.8            | 17.8        |
| LV end-diastolic volume       | mL           | 130.0         | 245.6      | 237.1       | 246.8           | 256.0       |
| RV end-diastolic pressure     | mmHg         | 4.0           | 3.3        | 16.8        | 4.8             | 20.5        |
| RV end-diastolic volume       | mL           | 121.0         | 106.4      | 252.8       | 125.6           | 277.5       |
| LV ejection fraction          | %            | 52            | 17         | 16          | 16              | 11          |
| RV ejection fraction          | %            | 56            | 39         | 15          | 31              | 10          |
| Pulmonary vascular resistance | Wood Units   | 2.0           | 2.0        | 3.1         | 9.0             | 9.0         |

\* The approximate calculation: pulmonary vascular resistance = (mean pulmonary arterial pressure – mean left atrial pressure)/cardiac output; LV: left ventricle; RV: right ventricle.



# Heart Assist Lab (BME Cornell)



## **References:**

1. Ursino M, Interaction between carotid baroregulation and the pulsating heart: a mathematical model, *Am J Physiol Heart Circ Physiol* 275: H1733-H1747, 1998
2. Santamore WP, Burkhoff D. Hemodynamic consequences of ventricular interaction as assessed by model analysis. *Am J Physiol* 1991; 260(1 Pt 2):H146-57.
3. Maughan WL, Sunagawa K, Sagawa K. Ventricular systolic interdependence: volume elastance model in isolated canine hearts. *Am J Physiol* 1987; 253(6 Pt 2):H1381-90.
4. Chung DC, Niranjana SC, Clark JW, Jr., et al. A dynamic model of ventricular interaction and pericardial influence. *Am J Physiol* 1997; 272(6 Pt 2):H2942-62.
5. Ferreira A., A rule-based controller based on suction detection for rotary blood pumps, University of Pittsburgh, 2007.